


NEW Sophi Products.



NEW Phaco Sets Designed for Higher Chamber Stability!

- Increased chamber stability¹
- Higher vacuum values applicable than with standard 20g/21g phaco set¹
- Recommended phaco set for 20g/21g venturi surgeons
- Suitable for peristaltic and venturi aspiration pump

Dynamic Phaco Set Sophi - 20G

Incision Size 2.5 – 2.75mm Aspiration Port 0.5mm

REF. 1401-3020 Single Use² REF. 1401-3120 Multi Use²

HighVac Phaco Set Sophi - 21G

Incision Size 2.2 – 2.4mm Aspiration Port 0.45mm

REF. 1401-2821 Single Use² REF. 1401-2921 Multi Use²



NEW Curved Phaco Sets!



Phaco Set Sophi - 20G Curved

Incision Size 2.5 – 2.75mm

REF. 1401-2720 Single Use² REF. 1401-3220 Multi Use²

Phaco Set Sophi - 21G Curved

Incision Size 2.2 - 2.4mm

REF. 1401-2721 Single Use² REF. 1401-3221 Multi Use²

Phaco Set Sophi - 23G Curved

Incision Size 1.8 - 2.0mm

REF. 1401-2723 Single Use² REF. 1401-3223 Multi Use²

NEW Sophi Eco!

- · Batteries optionally available
- · Batteries upgradeable at a later date
- Peristaltic aspiration

Sophi Eco

Active IOP Control Pump Peristaltic Pump Compressor Vitrectomy Battery-Ready Foot Pedal Sophi (1401-1300)

Optional:

License Battery Sophi (1401-3000) Tray Sophi (1401-0109) Wi-Fi Video Transmitter Sophi (1401-0199-300)

REF. 1401-0100-A



NEW Insole Foot Pedal Sophi!

- · Optimized for surgeons with shoe size 39 or smaller
- Insole retrofittable without assembling efforts



Insole Foot Pedal Sophi

Inner Length 246mm Inner Ankle Width 71mm Inner Depth 15mm

REF. 1401-9928-004



EC REP

This GmbH Millennium Park 9 6890 Lustenau, Austria

www.sophi.info



This AG
Widnauerstrasse 1
9435 Heerbrugg, Switzerland
T: +41 71 544 08 70
F: +41 71 544 10 86
info@thisag.ch
www.thisag.ch